

# PC-505B

## Patch Clamp

# patch clamp

*The lowest noise, switchable, resistive-feedback patch clamp amplifier currently available*



- Lowest noise — approaching theoretical limit
- Calibrated Cap Comp and Series R circuitry
- % Compensation circuitry
- Independent V hold and I hold controls
- Zap safety switch
- LED meter
- 3 year warranty

The PC-505B is the lowest noise, switchable, resistive-feedback patch clamp amplifier available. This model also has features of particular interest to those doing whole cell studies. The slow capacitance compensation circuitry has been combined into a single control and allows direct measurement of membrane capacitance. The companion Series R control displays the access resistance and the new % correction circuit compensates up to 90% of the access resistance. These and other features make the PC-505B an extremely capable amplifier.

### Switching Headstages

Two selectable feedback resistors in the headstages permit single channel and whole cell recording on the same cell. A 50 G $\Omega$  resistor is used in both switching models for low noise single channel recording with currents to 200 pA.

#### LC-201B Headstage (50 G $\Omega$ /500 M $\Omega$ )

This headstage with 500 M $\Omega$  feedback resistor will handle whole cell currents up to 20 nA.

#### HC-202B Headstage (50 G $\Omega$ /50 M $\Omega$ )

For larger whole cell currents, the 50 M $\Omega$  resistor in this headstage permits currents up to 200 nA.

### Bilayer Headstage

#### W4-205B Bilayer Headstage (50 G $\Omega$ /500 M $\Omega$ modified)

The 50 G $\Omega$  resistor headstage is modified for artificial bilayer capacitances up to 250 pF.

### Headstage Resistor Selection

Headstage feedback resistance is dynamically switched at the amplifier front panel. LED's indicate resistor selection and the corresponding multiplying factor applied to the current gain [Im] switch setting.

### Operating Modes

The PC-505B has three modes of operation: voltage clamp, zero current clamp, and current clamp.

#### V Clamp

In voltage clamp mode, the input range is  $\pm 1$  V. Active commands (V hold, junction and auto zero, test pulse, zap) and external inputs are scaled and summed at the headstage input. Capacity compensation and speed test are also active.

#### I<sub>o</sub>

Zero current is essentially a standby mode used to preset voltage hold or current hold levels before switching to voltage or current clamp. All commands are inactive with the exception of junction zero which, in this case, functions as an offset control for the electrode and tip potentials associated with the pipette.

#### I Clamp

Current clamp mode clamps the cell to a current level determined by the current hold setting and any external commands. Capacitance compensation is inactive in this mode.

### Commands

Voltage and current commands applied to the cell include: voltage and current hold, junction and auto zero, test pulse, speed test, and zap.

#### V & I Hold

Holding potentials and currents are set with separate controls eliminating the need to reset levels when switching between voltage and current clamp modes.

#### Junction & Auto Zero

Adjusting for offset potentials is performed with either the manual junction zero control or with auto zero. The offset potential is read on the meter or at the V<sub>c</sub> x10 output.

#### Test Pulse & Speed Test

Test Pulse and Speed Test are internally generated 50/60 Hz signals. Test pulse is attenuated by the command sensitivity and is useful for monitoring the formation of a gigaseal. The speed test signal is applied to the headstage input to allow for tuning the headstage response. A rear panel speed test switch allows for an external signal to be used.

## PC-505B

### Patch Clamp (continued)

#### Zap

Variable duration pulse used to rupture the cell membrane for whole cell recording.

Signals applied to the command input are attenuated at one of three levels with the command sensitivity selector.

#### Capacitance & Resistance Compensation: Fast Compensation

Stray capacitance between the input and electrode resistance is compensated with two pair of controls, C-Fast 1 and C-Fast 2. Amplitude and time constant of each pair is independently adjustable. Whole cell capacitance compensation is adjusted with the single control, C-Slow. Membrane capacitance is read from the C-Slow calibrated dial. The companion Series R control is used in conjunction with the C-Slow and its calibrated dial provides a reading of the access resistance. C-Slow may be disabled to view the uncompensated signal. The % Correction control is used to increase the command signal to compensate for the voltage drop across the access (Series R) resistance. Correction is adjustable up to 90%.

#### Outputs

Signals at the Im output are filtered at the selected cutoff frequency set with the 4-pole Bessel filter. Selecting Bypass presents the full bandwidth signal to the Im output. Filtered and unfiltered Im output is also available at the instrument rear panel. Additional outputs are the membrane voltage Vm x10 and the sum of the voltage commands Vc x10, both at x10 gain.

#### Voltammetry with PC-505B

PC-505B functions as an excellent low-noise potentiostat for voltammetric and other electrochemical measurements. In this mode, the V hold (electrode potential) is increased to a maximum of ±1 volt and the maximum external command signal to ±2 volts at electrode or sensor.

#### Specifications

##### Headstages:

###### LC-201B Headstage (50 GΩ/500 MΩ)

Single channel currents to 200 pA, whole cell currents to 20 nA.

###### HC-202B Headstage (50 GΩ/50 MΩ)

Single channel currents to 200 pA, whole cell currents to 200 nA.

###### W6-205B Bilayer Headstage (50 GΩ/500 MΩ modified)

For artificial bilayer capacitances up to 250 pF, currents to 20 nA.

##### Noise (referred to input)

Measured with an 8-pole Bessel filter, input open, 50 GΩ resistor:

|             |              |
|-------------|--------------|
| DC to 1 kHz | 0.035 pA RMS |
| DC to 5 kHz | 0.150 pA RMS |

|           |        |
|-----------|--------|
| Bandwidth | 25 kHz |
|-----------|--------|

##### Voltage Clamp Commands:

|                     |   |
|---------------------|---|
| Command In BNC      | ±10 V Max, AC or DC, applied to input                                   |
| Voltage Hold        | ±200 mV Max with 10-turn control  |
| Junction Zero       | ±100 mV Max with 10-turn control  |
| Internal Test Pulse | 1 V 50/60 Hz (line freq.) square wave attenuated by Command Sensitivity |
| Command Sensitivity | x0.1, x0.01, and x0.001   |
| Zap                 | 1.0 V Pulse, adjustable duration from 0.1 to 10 msec                    |

##### Current Clamp Commands:

|                     |  |
|---------------------|--|
| Command In          | ±1000 pA max with Command Sensitivity @ x0.1<br>±100 pA max with Command Sensitivity @ x0.01<br>±10 pA max with Command Sensitivity @ x0.001 |
| Current Hold        | ±1 nA with 10-turn control   |
| Internal Test Pulse | 1 nA 100 Hz square wave through Command Sensitivity  |
| Command Sensitivity | x 0.1, x 0.01, and x 0.001   |

#### Specifications (continued)

##### Fast Capacitance Compensation (Voltage Mode):

|          |                              |
|----------|------------------------------|
| C Fast 1 | 0.1 to 1.75 μsec, 0 to 5 pF  |
| C Fast 2 | 0.33 to 8.5 μsec, 0 to 15 pF |

##### Whole Cell Compensation:

|              |                               |
|--------------|-------------------------------|
| C Slow       | 0-100 pF with 10 turn control |
| Series R     | 0-10 MΩ with 10 turn control  |
| % Correction | 0-90% of Series R             |

|                  |   |
|------------------|---|
| Leak Subtraction | 50 GΩ Headstage Resistor ∞ to 50 GΩ<br>500 MΩ Headstage Resistor ∞ to 500 MΩ<br>50 MΩ Headstage Resistor ∞ to 50 MΩ |
|------------------|---|

Front Panel Outputs: Im (membrane current), selected in the range of:  
Gains of 0.05 to 10 mV/pA with 50 MΩ headstage resistor  
Gains of 0.5 to 100 mV/pA with 500 MΩ headstage resistor  
Gains of 5 to 1000 mV/pA with 50 GΩ headstage resistor

|                    |   |
|--------------------|---|
| Vc x10             | Summation of all commands amplified by 10                           |
| Vm x10             | Membrane voltage amplified by 10                                    |
| Im Low-Pass Filter | 0.05, 0.1, 0.2, 0.5, 1, 2, 5, 10, 20 kHz and Bypass (4-Pole Bessel) |

##### Panel Meter:

|               |   |
|---------------|---|
| Junction Zero | ±199.9 mV                                       |
| Vc + h In     | Sum of all commands and V Hold, ±199.9 mV       |
| Vc            | Sum of all commands, ±199.9 mV                  |
| Vm            | Membrane voltage (current clamp mode) ±199.9 mV |
| Im            | Membrane current, ±1999 pA                      |
| RMS Noise     | 1.999 pA  |

##### Rear Panel Outputs:

|                   |   |
|-------------------|---|
| Gain Telegraph    | From 0.5 to 7.0 V in 0.5 V steps.*                        |
| Filter Telegraphs | From 0.2 to 2.0 V in 0.2 V steps.*                        |
| Im/Vm Telegraph   | Logic levels, V-Clamp=1, I-Clamp=0                        |
| Sync Out          | Signal for synchronizing an oscilloscope to internal test |

##### Power Requirements

110 to 130 or 220 to 250 VAC, 50/60 Hz, 15 VA

##### Physical Dimensions:

|              |   |
|--------------|---|
| Main Unit    | 8.9 x 43.2 x 30.5 cm, H x W x D                 |
| Headstage    | 1.9 x 3.5 x 5.7 cm, H x W x L, with 1.8 m cable |
| Mounting Rod | 6.3 mm D x 6.3 cm cm L                          |

Shipping Weight 11.4 kg

Warranty Three years, parts and labor

\* Compatible with Axon pClamp and Heka PatchMaster.

#### Order # Model Product

##### Line operating voltage if other than 100-130 VAC.

|            |           |   |
|------------|-----------|---|
| W4 64-0000 | PC-505BLC | Patch Clamp PC-505B with LC-201B Headstage*         |
| W4 64-0001 | PC-505BHC | Patch Clamp PC-505B with HC-202B Headstage*         |
| W4 64-0002 | PC-505BHB | Patch Clamp PC-505B with HB-205B Bilayer Headstage* |

\* Supplied with model cell and rack mount hardware.

#### Additional/Replacement Headstages

|            |         |                             |
|------------|---------|-----------------------------|
| W4 64-0004 | LC-201B | 50 GΩ/500 MΩ Headstage      |
| W4 64-0005 | HC-202B | 50 GΩ/50 MΩ Headstage       |
| W4 64-0006 | HB-205B | 50 GΩ Headstage for Bilayer |

#### Electrode Holders for PC-501A and PC-505B

|            |          |                                 |
|------------|----------|---------------------------------|
| W4 64-0821 | QSW-A10P | Straight Holder 1.0 mm glass OD |
| W4 64-0822 | QSW-A12P | Straight Holder 1.2 mm glass OD |
| W4 64-0823 | QSW-A15P | Straight Holder 1.5 mm glass OD |
| W4 64-0978 | QSW-A17P | Straight Holder 1.7 mm glass OD |
| W4 64-0824 | QSW-A20P | Straight Holder 2.0 mm glass OD |