


other coolers/warmers


PHC Series

Heater/Cooler Jackets

Maintains temperature in both perfused and static baths



Model PHC-3 shown with mounting platform and mounted in a Nikon stage adapter




Individual heater/cooler jackets shown without mounting platform

- Designed for Series 20 chambers
- Optimized for the SC-20 In-line Solution Heater/Cooler
- Jackets available for both upright and inverted microscopes
- Includes Series 20 platform

The Warner PHC Heater/Cooler Jackets are designed to bring heating and cooling to our classic Series 20 Imaging and Recording Chambers. Accurate temperature control from 5° to 50°C can be achieved using the PHC jackets in concert with the SC-20 Dual In-line Solution Heater/Cooler. Heated or chilled water flows from the SC-20 into a PHC jacket which is in direct thermal contact with the chamber bottom coverslip.

The PHC-1 is used for upright microscopes and provides a thermal barrier between the chamber-forming coverslip and the local environment. The PHC-2 and PHC-3 are designed for inverted microscopes and provide either rectangular or round openings.

Heater/Cooler Jackets are provided with a mounting platform, which replaces the standard platforms used with Series 20 chambers. The platform functions as a base for the jacket/chamber sandwich and provides the clamping pressure to make a tight seal. Mounting platforms are machined from black Delrin and are compatible with all Series 20 stage adapters, see pages 64 to 69.



Assembly of Series 20 chamber with cooling jacket and stage adapter

Model	Aperture Size	For Chamber Model
PHC-1	17.0 x 37 mm	RC-22/22C/24N/26/26G/26GLP/26Z RC-27/27N/27NE/28/RC-27L/RC-29
PHC-2	15 mm diameter	RC-22/22C/24N/26/26G/26Z/26GLP
PHC-3	8.0 x 25 mm	RC-27/27N/27NE

Order #	Model	Product
W4 64-0354	PHC-1	Heater/Cooler Jacket, Upright
W4 64-0355	PHC-2	Heater/Cooler Jacket, Inverted
W4 64-0356	PHC-3	Heater/Cooler Jacket, Inverted

Temperature Control

other coolers/warmers