

# microelectrode holders

## MP Series

### MP Series Electrode Holders

Microelectrode holders for every application and manufacturer



- Microinjection

#### Pressure Ports

A	Accepts 1.5 mm (1/16") tubing
B	Accepts 2.4 mm (3/32") tubing*
ML	Accepts female Luer fitting
FL	Accepts male Luer fitting
T*	Accepts 1/4"-28 male thread fitting

#### Specifications

Body Material	Acrylic
Body Style	Straight
Port	1/16" barbed (1/16" tubing ID), 3/32" barbed (3/32" tubing ID), male Luer, female Luer and 1/4"-28 female threaded
Glass Sizes	1.0, 1.2, 1.5 and 2.0 mm
Handle	6.3 mm diameter x 6.3 cm long, threaded to attach to holder.

#### MP SERIES Holders

Order #	Body Style	Port	Glass OD (mm)	Model	
W4 64-1261	Straight	1/16" barb	1.0	MP-S10A	
W4 64-1262			1.2	MP-S12A	
W4 64-1263			1.5	MP-S15A	
W4 64-1264			2.0	MP-S20A	
W4 64-1265		3/32" barb*	1.0	MP-S10B	
W4 64-1266			1.2	MP-S12B	
W4 64-1267			1.5	MP-S15B	
W4 64-1268			2.0	MP-S20B	
W4 64-1273		male Luer	1.0	MP-S10ML	
W4 64-1274			1.2	MP-S12ML	
W4 64-1275			1.5	MP-S15ML	
W4 64-1276			2.0	MP-S20ML	
W4 64-1269		female Luer	1.0	MP-S10FL	
W4 64-1270			1.2	MP-S12FL	
W4 64-1271			1.5	MP-S15FL	
W4 64-1272			2.0	MP-S20FL	
W4 64-1278		threaded	1.0	MP-S10T	
W4 64-1279			1.2	MP-S12T	
W4 64-1280			1.5	MP-S15T	
W4 64-1281			2.0	MP-S20T	

\* Shaft Diameter

\*\* Threads/Inch