

# Tubing and Connector Kit, Syringe Needles

## Parts and accessories




Warner Instruments offers a complete tubing and connector kit. This kit is ideal for use with Warner's Imaging and Recording Chambers. The kit is comprised of barbed and Luer fittings, blunt end needles, and an assortment of PE and C-Flex tubing. All barbed and Luer fittings are made from polypropylene and include tube to tube, reducing, Y-, T-, Luer-to-Luer, and Luer-to-barb adapters.

Comes in a convenient plastic storage box.

Order #	Model	Product
W4 64-1565	KIT-1	Fitting and Tubing Kit

Components listed at right are included in the tubing and connector kit. Components also sold separately.



Order #	Product
W4 64-1566	Tube Fitting Barb 1/8" to 1/16", pkg. of 10
W4 64-1567	Tube Fitting Barb 1/8" to 1/16", pkg. of 10
W4 64-1568	Tube Fitting Barb 1/8" to 1/8" in, pkg. of 10
W4 64-1569	Tube Fitting Tee Barb 1/16", pkg. of 10
W4 64-1570	Tube Fitting Tee Barb 1/8", pkg. of 10
W4 64-1571	Tube Fitting Y Barb 1/16", pkg. of 10
W4 64-1572	Tube Fitting Y Barb 1/8", pkg. of 10
W4 64-1573	Tube Fitting Barb 1/16" to Luer Male, pkg. of 10
W4 64-1574	Tube Fitting Barb 1/8" to Luer Male, pkg. of 10
W4 64-1575	Tube Fitting Barb 1/8" to Luer Female, pkg. of 10
W4 64-1576	Tube Fitting Barb 1/8" to Luer Female, pkg. of 10
W4 64-1577	Tube Fitting Luer Male to Luer Female, pkg. of 10
W4 64-1578	Tube Fitting Luer Tee Female, pkg. of 10
W4 64-1579	Tube Fitting Luer Male to Luer Male, pkg. of 10
W4 64-1580	Tube Fitting Luer Female to Luer Female, pkg. of 10
W4 64-1581	Tube Fitting Luer Female Plug, pkg. of 10
W4 64-1582	Tube Fitting Luer Male Plug, pkg. of 10
W4 64-0141	LPE-50, Luer to PE-50 tubing adapter, pkg. of 8
W4 64-0166	CFL-6, C-Flex tubing (1/32" ID x 6 ft) 3/32 OD mm, 6 ft

### Blunt Needle Plastic Hub

W4 64-1489	SN-18, 18G, 0.5", pkg. of 12
W4 64-1490	SN-23, 23G, 0.5", pkg. of 12
W4 64-1583	SN-30, 30G, 0.5", pkg. of 12

### Polyethylene Tubing

W4 64-0750	PE-10/10, 0.28 ID x 0.61 OD mm, 10 ft
W4 64-0752	PE-50/10, 0.58 ID x 0.97 OD mm, 10 ft
W4 64-0754	PE-90/10, 0.86 ID x 1.27 OD mm, 10 ft
W4 64-0755	PE-160/10, 1.14 ID x 1.57 OD mm, 10 ft